# Ivotion Base

<b>[en]</b> <b>Instructions for Use</b> Denture base material for the CAD/CAM technique	<b>[sl]</b> <b>Navodila za uporabo</b> Material osnove zobne proteze za tehniko CAD/CAM
<b>[de]</b> <b>Gebrauchsinformation</b> Prothesenbasismaterial für die CAD/CAM-Technik	<b>[hr]</b> <b>Upute za uporabu</b> Materijal za bazu proteze za CAD/CAM tehniku
<b>[fr]</b> <b>Mode d'emploi</b> Matériau pour base prothétique pour la technique CAD/CAM	<b>[cs]</b> <b>Návod k použití</b> Materiál na báze snímatelných náhrad pro CAD/CAM techniku
<b>[it]</b> <b>Istruzioni d'uso</b> Materiale per basi protesiche per la tecnica CAD/CAM	<b>[sk]</b> <b>Návod na používanie</b> Základný materiál na zubné protézy pre techniku CAD/CAM
<b>[es]</b> <b>Instrucciones de uso</b> Material para bases de dentadura mediante la técnica CAD/CAM	<b>[hu]</b> <b>Használati utasítás</b> Protézisalapanyag CAD/CAM technikához
<b>[pt]</b> <b>Instruções de Uso</b> Material para base de próteses para técnica CAD/CAM	<b>[sr]</b> <b>Упутство за употребу</b> Материјал за основу протетске надокнаде за CAD/CAM технику
<b>[sv]</b> <b>Bruksanvisning</b> Protesbasmaterial för CAD/CAM-tekniken	<b>[mk]</b> <b>Упатство за употреба</b> Материјал за база на протезите за CAD/CAM техника
<b>[da]</b> <b>Brugsanvisning</b> Protesebasismateriale til CAD/CAM-teknikken	<b>[bg]</b> <b>Инструкции за употреба</b> Материал за основа на протези за CAD/CAM техника
<b>[fi]</b> <b>Käyttöohjeet</b> Hammasproteesin runkomateriaali CAD/CAM-menetelmään	<b>[sq]</b> <b>Udhëzimet e përdorimit</b> Material baze proteze për teknikën CAD/CAM
<b>[no]</b> <b>Bruksanvisning</b> Protesebasismateriale til CAD/CAM-teknologi	<b>[ro]</b> <b>Instrucțiuni de utilizare</b> Material de bază de proteză pentru tehnica CAD/CAM
<b>[nl]</b> <b>Gebruiksaanwijzing</b> Basismateriaal voor gebitsprothesen voor de CAD/CAM-techniek	<b>[uk]</b> <b>Інструкція з використання</b> Базовий матеріал протезів для систем CAD/CAM
<b>[el]</b> <b>Οδηγίες Χρήσης</b> Υλικό βάσης οδοντοστοιχίας για την τεχνική CAD/CAM	<b>[et]</b> <b>Kasutusjuhend</b> Hambaproteesi alusmaterjal CAD-/CAM-tehnika jaoks
<b>[tr]</b> <b>Kullanım Talimatları</b> CAD/CAM tekniği için protez kaide materyali	<b>[lv]</b> <b>Lietošanas instrukcija</b> Zobu protēžu bāzes materiāls CAD/CAM metodei
<b>[ru]</b> <b>Инструкция по применению</b> Материал для базисов съёмного протеза при использовании технологии CAD/CAM	<b>[lt]</b> <b>Naudojimo instrukcija</b> Protezų pagrindo medžiaga, skirta CAD / CAM technikai
<b>[pl]</b> <b>Instrukcja stosowania</b> Materiał na płyty protez do techniki CAD/CAM	

CE 0123

Rx ONLY

Date information prepared:  
2024-01-23 / Rev. 0

**Manufacturer:**  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
[www.ivoclar.com](http://www.ivoclar.com)

**ivoclar**

## 1 Intended use

### Intended purpose

Fabrication of bases for removable dentures

### Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

### Intended users

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- Denturists (fabrication of restorations in the dental laboratory)

### Special training

No special training required.

### Use

For dental use only.

### Description

Ivotion Base are impact-resistant PMMA discs for the CAD/CAM technique to fabricate denture bases in removable denture prosthetics using a subtractive manufacturing process.

The discs are made of high-impact modified PMMA material.

They are available in a variety of shades to faithfully reproduce the individual characteristics of the patient.

### Technical data

	According to ISO 20795-1:2013*	Typical materials properties
Flexural strength	≥ 65 MPa	69 MPa
Flexural modulus	≥ 2000 MPa	2096 MPa
Fracture toughness K <sub>max</sub>	≥ 1.9 MPam <sup>1/2</sup>	2.1 MPam <sup>1/2</sup>
Overall work of fracture W <sub>f</sub>	≥ 900 J/m <sup>2</sup>	1224 J/m <sup>2</sup>
Water absorption	≤ 32 µg/mm <sup>3</sup>	18.3 µg/mm <sup>3</sup>
Water solubility	≤ 1.6 µg/mm <sup>3</sup>	0.5 µg/mm <sup>3</sup>
Residual MMA content	≤ 4.5%	2.5 %

\* Currently, there is no standard that defines the physical requirements for machinable denture base materials processed with milling systems. The product complies with the requirements of ISO20795-1:2013 for heat-curing denture base resins. The residual MMA content complies with the requirements of self-curing denture base materials.

### Indications

- Partial edentulism in the anterior and posterior region
- Complete edentulism

### Types of restorations:

- Complete dentures
- Partial dentures
- Implant-supported dentures
- Tooth-supported dentures

### Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

### Limitations of use

- Observe the required minimum wall thickness (2 mm). Failure to observe the minimum thickness may lead to failure (e.g. fracture of the restoration).
- Only mill the discs in a compatible CAD/CAM system.
- Staining / layering with SR Nexco (light-curing, MMA-based conditioner for the bonding of light-curing materials to PMMA materials (heat-curing or self-curing polymers and denture teeth)).

- Avoid overheating of the material when finishing the objects.
- Avoid contact with solvents.
- Protect the materials from direct sunlight.
- Failure to observe the stipulated limitations of use may lead to failure.
- The final denture is for use by one patient only.

### Side effects

In individual cases, allergic reactions to methyl methacrylate materials have been reported.

### Interactions

There are no known interactions to date.

### Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

### Composition

PMMA, high-impact modified PMMA, pigments

## 2 Application

**NOTICE! Avoid overheating the material through improper use of the milling tools, polishing brushes and steam cleaners.**

### Designing the denture

**NOTICE! Observe the following specifications:**

- Wall thickness for denture base: min 2 mm.
  1. Design the denture.

### Fabricating the denture base

**NOTICE! Use a compatible milling system with a standard holder of 98.5 mm.**

**NOTICE! Only use milling tools recommended for the particular milling system.**

**NOTICE! Make sure that the disc holder is securely fixed.**

1. Mill the denture base.

### Preparing the denture base

1. Detach the holding bars. Use a suitable tungsten carbide bur.
2. Smooth out the attachment points.

### Optional: Finishing the denture base

**NOTICE! Pay special attention to the margins of the tooth cavities, the interdental spaces and the anatomical surface structure.**

1. Finish the denture base. Use a cross-cut or diamond tungsten carbide bur.

### Checking the fit

1. Check the fit of the teeth in the denture base.
2. In case of interfering contacts: Adjust the tooth. Use a cross-cut or diamond tungsten carbide burs.

### Bonding the teeth to the denture base

1. Bond the teeth on the denture base. Use a suitable bonding system (e.g. Ivotion Bond).  For guidance regarding processing, consult the Instruction for Use of the product.

### Finishing/polishing the denture

1. Finish the denture. Use cross-cut or diamond tungsten-carbide burs.

**NOTICE! Gloss varnish must not be used.**

2. Optional: Characterize the denture. Use the lab composite SR Nexco.  For guidance regarding processing, consult the Instruction for Use of the product.
3. Polish the denture. Use a goat-hair brush and pumice.
4. Polish the denture to a high gloss. Use a cotton or leather polishing buff and polishing paste.

### Cleaning the denture

1. Clean the denture under running water. Use a tooth brush.

### Repairing/extending/relining the denture

1. Repair/extend/reline the denture. Use a conventional self-curing denture resin (e.g. ProBase® Cold)  For guidance regarding processing, consult the Instruction for Use of the product.

### 3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, [www.ivoclar.com](http://www.ivoclar.com) and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: [www.ivoclar.com/elfu](http://www.ivoclar.com/elfu)
- Basic UDI-DI: 76152082APROS007W7

### Warnings

- Observe the Safety Data Sheet (SDS) (available at Ivoclar Vivadent AG [www.ivoclar.com](http://www.ivoclar.com)).
- Do not inhale grinding dust. Use an extraction unit and wear a face mask.
- Wear protective goggles during finishing.
- Processing the material may result in sharp edges. There is a risk of injury.

### Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Ingestion of fragments
- Occurrence of denture stomatitis

### 4 Shelf life and storage

- Protect the materials from direct sunlight.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

### 5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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- Perskaitykite saugos duomenų lapą (SDL) (jį rasite Ivoclar Vivadent AG [www.ivoclar.com](http://www.ivoclar.com)).